Benchmark Reference Standards

PD-L1 Quantitative Control

Cell line controls for immunohistochemistry and in situ hybridization.

Research Use Only

CC-BRS2023 2 Slide Pack CC-BRS2024 5 Slide Pack

CC-BRS2025 X1 Cell microarray block

Application

This product is suitable for use in immunohistochemistry and in situ hybridization.

Materials

Four formalin fixed paraffin embedded cell lines with a quantitative range of expression for PD-L1.


Cell line A: Breast ductal carcinoma

Cell line B: Osteosarcoma Cell line C: Fibrosarcoma

Cell line D: T-cell non-Hodgkin lymphoma

Cells are fixed in 10% neutral buffered formalin and paraffin wax embedded.

Sections are cut at 4 μ m, mounted on positively charged slides and baked overnight at 37°C. Cell microarrays (CMA) contain cores that are 1.5-2mm in diameter and 3-3.5mm in length. It is possible to obtain up to 400 sections depending on thickness and batch cutting.





Expression Profile

Cell Line	PD-L1
А	Negative
В	Low
С	Medium
D	High

Storage and Handling

Store at 2-8°C. Do not freeze (for expiration date please see the product label)

Warnings and Precautions

- 1. The product is intended for research use only.
- It is the responsibility of the end user to determine suitability with their reagents and procedures within their laboratory.
- 3. Do not use after expiration date printed on product labels. The user must validate any storage conditions other than those specified in the package insert.

Troubleshooting

For further help please feel free to contact cancercontrols.com at service@cancercontrols.com.



For updates and additional product information please visit: cancercontrols.com